

leon *mri*



The **LEON mri** has been specially developed for use in the MRI area with field strengths of up to 40 millitesla. It can be operated on 1.5 tesla und 3.0 tesla MRI systems. Different coloured alarms and messages indicating the priority can be easily recognised due to LED lights, which are visible from all sides, also outside the MR field.

Power spectrum of the **LEON mri**:

- Both the power spectrum and the operation of the **LEON mri**, match exactly, without any restrictions, those of the already well-established **LEON**.
- A further major logistic advantage can be found in the compatibility of the compact circuit system with the **LEON plus**, **LEON** and **LEON mri** systems.
- The position to MRI is monitored by an integrated magnetic field strength monitor.

