

# PEARL SpO<sub>2</sub>

High quality probes for  
Medlab Pulse Oximeters



Medlab PEARL SpO<sub>2</sub> probes can be used with all Medlab pulse oximeters currently sold. They are the result of more than 15 years of experience in developing SpO<sub>2</sub> sensors.

By using only the highest quality cabling, a Kevlar re-enforcement and a specially designed cable strain relief, we are able to offer all our reusable probes with a two year warranty.

The standard cable length of the sensors is 140 cm, extension cables are optionally available.

Please see the table on the right side for more information. We deliver a full line of probes for different applications.

By using the PEARL probes, your budget for pulse oximetry is easy to plan: You can be sure that for 24 months, there will be no additional costs because of broken probes\*.



Medlab medizinische Diagnosegeräte GmbH  
Helmholtzstrasse 1, 76297 Stutensee near Karlsruhe, Germany  
Tel. +49 (0) 7244 74110-0, FAX +49 (0) 7244 74110-28  
sales@medlab.eu, www.medlab.eu

**medlab**  
medizinische Diagnosegeräte GmbH

## Ordering Info

### Reusable Probe Selection:

Type	Patient Weight	Order Number
Fingerclip P-200	> 20 kg	00535
Fingerclip PS-200	> 20 kg	00537
Oxiflex large R-200	> 20 kg	00550
Oxiflex small RS-200	> 10 kg	00570
Wrap Probe WR-200	> 1 kg	00532
Y Probe Y-200	> 1 kg	00526
Tongue probe V-200	> 20 kg	00440

### Disposable Probe Selection:

Type	Patient Weight	Order Number
Adult, Microfoam™	> 20 kg	00520
Adult, Plaster	> 20 kg	00540
Child, Microfoam™	8 kg - 30 kg	00521
Child, Plaster	8 kg - 30 kg	00541
Infant, Microfoam™	5 kg - 15 kg	00522
Infant, Plaster	5 kg - 15 kg	00542
Neonate, Microfoam™	< 5 kg	00523
Neonate, Plaster	< 5 kg	00543



P-200  
Fingerclip probe -  
standard delivery, 140 cm



PS-200  
Small Fingerclip probe, 140 cm



R-200  
Oxiflex large finger probe, 140 cm



RS-200  
Oxiflex small finger probe,  
140 cm



WR-200  
Wrap probe, 140 cm



Y-200  
Universal Y-probe, 120 cm



V-200  
Tongue probe - for veterinarian  
applications, 140 cm



Disposable probes  
Delivery:  
Package of 24 pieces

\*Medlab PEARL probes can only be used for current Medlab pulse oximeters. If in doubt, contact your dealer or Medlab directly. Mechanical damages of the probe, as, for example, cutting the cable are exempted from warranty. PEARL® is a registered trademark of Medlab GmbH.