

1. TOF-Watch devices and accessories:

| Content: | | |
|---|---|-------------|
| Organon article no. | Product description | Page |
| Devices: | | |
| 79950161 | TOF Watch sales package | 2 |
| 79950160 | TOF Watch S sales package | 2 |
| 79950162 | TOF Watch SX sales package | 2 |
| 79950157 | TOF Watch SX + Interface (Philips Medical) | 3 |
| General accessories: | | |
| 79950074 | Negative stimulation lead (5 pcs.) | 4 |
| 79950075 | Positive stimulation lead (5 pcs.) | 4 |
| 79950163 | Acc. transducer lead (TOF-Watch / TOF-Guard) | 4 |
| 79950215 | Main TOF Watch cable | 5 |
| 79950006 | Local anesthesia cable for TOF-Watch | 5 |
| 79950152 | Thumb adapter acceleration transducer - 50 pcs. | 5 |
| 79950207 | Eye adapter acceleration transducer - 50 pcs. | 6 |
| 79950014 | TOF-Watch IV-Pole clamp | 6 |
| 79950216 | Battery lid for TOF-Watch (25 pcs.) | 6 |
| 79950209 | Hand adapter for acc. transducer | 7 |
| 79950217 | Strap for Hand adapter (10 pcs.) | 7 |
| 79950039 | Electrodes Neotrode (3 x 100 pcs.) | 7 |
| Accessories for TOF-Watch / TOF-Watch S: | | |
| 79950214 | Complete TOF-Watch cable | 8 |
| Accessories for TOF-Watch SX | | |
| 79950213 | Complete TOF-Watch SX cable | 9 |
| 79950151 | Thermistor lead (5 pcs.) | 9 |
| TOF-Link and accessories: | | |
| 79950155 | TOF-Link USB sales package | 10 |
| 79950154 | TOF Link cable 3 m | 10 |
| 79950208 | TOF-Link cable 0.25m | 10 |
| 79950219 | TOF-Watch SX Interface RS232 adapter | 11 |

a. Devices:

ART. NO.
DESCRIPTION

PICTURE

79950161

TOF-Watch

Containing:

1 Instructions for use.

1 TOF-Watch.

1 TOF-Watch Patient cable
(79950214)



79950160

TOF-Watch S

Containing:

1 Instructions for use.

1 TOF-Watch S

1 TOF-Watch Patient cable
(79950214)



79950162

TOF-Watch SX

Containing:

1 Instructions for use

1 TOF-Watch SX

1 TOF-Watch SX Patient cable
(79950213)



ART. NO.
DESCRIPTION

PICTURE

79950157
TOF-Watch SX + interface (Philips)

Containing:

- 1 Instructions for use
- 1 TOF-Watch SX
- 1 TOF-Watch SX Patient cable (79950213)
- 1 TOF-LINK RS 232 interface
- 1 TOF-LINK interface cable (0.25 m) (79950208)



b. General accessories:

ART. NO.
DESCRIPTION

PICTURE

79950074
Negative stimulation lead (5 pcs.)

Part of all TOF-Watch sales packages and TOF-Watch complete cables.



79950075
Positive stimulation lead (5 pcs.)

Part of all TOF-Watch sales packages and TOF-Watch complete cables.



79950163
TOF-Watch SX / TOF Guard acceleration transducer

Part of all TOF-Watch sales packages and TOF-Watch complete cables.



ART. NO.
DESCRIPTION

PICTURE

79950215

TOF-Watch main cable

Part of all TOF-Watch sales packages and TOF-Watch complete cables.



79950006

Stimulation cable (1.5 m) for needle electrode (LA)

One lead to be connected to a surface electrode and the 2 mm (diameter) banana to be connected to a Local Anesthesia needle.



79950152

Thumb adapter acceleration transducer - 50 pcs.

For attaching the acceleration transducer to the thumb (without plaster). Using this adapter avoids damage to the cable when detaching the cable.



ART. NO.
DESCRIPTION

PICTURE

79950207

**Eye adapter acceleration
transducer - 50 pcs.**

For attaching the acceleration
transducer to the orbicularis oculi or
corrugator supercilii.
Box of 50 pcs. eye adapter includes
50 pcs ready to use double sided
adhesives.



79950014

**Mounting bracket (clamp for IV-
pole)**

Fits all TOF-Watch devices



79950216

Battery lid for TOF-Watch (25 pcs)



ART. NO.
DESCRIPTION

PICTURE

79950209
Handadapter for acceleration transducer.

Containing:

- 1 Handadapter
- 2 Straps for handadapter (79950014)



79950217
Strap for Handadapter (10 pcs)


Part of article number 79950209
(2 pcs.)



79950039
Electrodes Neotrode (3 x 100 pcs.)



c. Accessories for TOF-Watch / TOF-Watch S:

| ART. NO. DESCRIPTION | PICTURE |
|--|--|
| <p>79950214 TOF-Watch complete cable.</p> <p>Part of article number 79950161 and 79950162. <u>Containing:</u> 1 TOF-Watch main cable (79950215) 1 Negative stimulation lead (79950074) 1 Positive stimulation lead (79950075) 1 acceleration transducer (79950163)</p> |  |

d. Accessories for TOF-Watch SX:

ART. NO.
DESCRIPTION

PICTURE

79950213
TOF-Watch SX complete cable.

Part of article number 79950162 and
79950157


Containing:

- 1 TOF-Watch main cable
(79950215)
- 1 Negative stimulation lead
(79950074)
- 1 Positive stimulation lead
(79950075)
- 1 acceleration transducer
(79950163)
- 1 TOF-Watch SX thermistor
(79950151)



79950151
TOF-Watch SX thermistor (5 pcs).

Part of article number 79950162 and
79950157 (1. pcs.)



e. TOF-Link and accessories:

ART. NO.
DESCRIPTION

PICTURE

79950155
TOF-LINK USB interface incl.
TOF-MONITOR program

Optical to USB adapter, cable and TOF-Watch SX Monitor program facilitates online TOF-Watch SX data capture on a Computer (desktop or laptop)

Containing:

- 1 CD with TOF-Watch SX Monitor program.
- 1 TOF-LINK USB interface
- 1 TOF-LINK interface cable (3m) (79950154)



79950154
TOF-LINK interface cable (3 m)

Part of article number 79950155



79950208
TOF-LINK interface cable (0,25 m)

Part of article number 79950157



ART. NO.
DESCRIPTION

PICTURE

79950219
TOF-Watch SX Interface RS232
adapter

Part of article number 79950157

